



Mopec Portable Bone Dust Vacuum



[MOPEC.COM/PORTABLE-BONE-DUST-VACUUM](https://www.mopec.com/portable-bone-dust-vacuum)

Portable Bone Dust Vacuum - BD050

Designed to fit the Mopec 810 and Mopec 1000 Autopsy Saws, the Mopec Portable Bone Dust Vacuum features a lightweight construction and convenient handle that provides easy mobility for the user. The self contained unit conveniently stores replacement autopsy saw blades, keeping everything you need in one place.



SPECIFICATIONS

Construction

- Dust collection unit with (1) pre-installed HEPA filter and grounded receptacle for an autopsy saw
- (1) 6' vacuum hose with end fittings for the vacuum inlet and vacuum nozzle
- (1) vacuum nozzle for use with autopsy saw
- (2) storage caps for the filter cartridge
- (6) hose straps
- (1) 10' power cord

Controls

- When autopsy saw is connected to the vacuum receptacle, the saw will only operate when the vacuum switch is in the "on" position.

Dimensions

- 7.5" W x 9.75" H x 19.5" L
- Weight: 10 lbs.

Utilities / Power

- 115V / 50-60 Hz
- 6.75 amps

Performance

- Disposable HEPA filter cartridge is easy to replace and clean
- Self-contained unit conveniently stores replacement autopsy saw blades
- Lightweight construction and convenient handle provides easily mobility

Compatibility

- Compatible with (BD810) Mopec 810 Autopsy Saw and (BD040) Mopec 1000 Autopsy Saw

Options & Accessories

- BD051: Replaceable HEPA Filter
- BD020: Saw Mount Bracket Kit
- BD021: Saw Attachment Y Bracket
- BD022: BD810 Vacuum Adapter



www.mopec.com
800.362.8491
800 Tech Row
Madison Heights, MI 48071