

Mopec *Rhythm*
MANUAL MICROTOME



DESIGNED WITH RELIABILITY IN MIND.

Mopec Rhythm Manual Microtome

Built specifically for use in histology and pathology laboratories, the Mopec Rhythm features a manual operation system and durable construction, giving the user precise control for the sectioning process. Designed with reliability in mind, this simple-to-use microtome is built to be more efficient and consistent in microtomy.

CONSISTENT BLOCK ALIGNMENT

Quick Cassette Clamp

Easily change out blocks with the easy-to-use, large lever that quickly releases a block for sectioning. The built-in level ensures consistent block alignment every time, providing the user with a quick reference.



ENHANCED USER CONTROL


Ergonomic Hand Wheel

Control the vertical movements of the blade with the ergonomic hand wheel. Built for enhanced user safety and comfort, the hand wheel allows the user to accurately adjust the thickness of the sections being cut, producing high-quality, uniform sections for analysis.

LONGER BLADE LIFE

2-in-1 Blade Holder

Easily switch between tissue types, blade supply, or technician preferences with the 2-in-1 blade holder. The simple click-into-place positioning system allows the user to easily slide the blade holder into to the left or right, allowing you to utilize both sides of the blade. The German-imported Kipp levers ensure a secure locking mechanism of the blade holder, giving you quality sections every time.



CLEAN SECTIONING

Manual Coarse Feed

The manual coarse feed allows for a more natural control over the specimen head.

Waste Tray

Easily catch sectioning waste with the lightweight, simple, and sleek waste tray. Simply slide the tray under the front of the microtome and ensure an easy cleanup.