

# Mopec Swordfish Autopsy Saw

Designed with the user in mind, the Mopec Swordfish Submersible Autopsy Saw is equipped with state-of-the-art risk prevention features. It boasts advanced technology and a sealed housing specifically constructed for autopsy use. This cutting-edge tool is available exclusively from Mopec.



## SPECIFICATIONS

### Construction

- Includes Saw Handpiece, Power Supply, and (2) Power Cords
- Completely submersible for up to 15 minutes to ensure effective cleaning
- Hermetically sealed to prevent shock from liquid contact
- Single press on/off button for use
- Device is double pole fused for added protection

### Weight & Dimensions

- Saw without power cord: 2.8 lbs/1.3 kg
- Power Supply: 3.7 lbs/1.7kg
- Saw: 9.5" L x 2.5" W
- Power Supply: 10" L x 4.75" W x 4.5" H
- Saw Power Cord: 12.5 ft. L
- Power Supply Cord: 8.5 ft. L

### Materials

- Saw Handpiece: Anodized aluminum body with 5 bearing sealed mechanism

### Utilities / Power

- US: 115 V
- International: 230 V

### Performance

- Oscillation: 13,200 CPM
- Operates on 40V to ensure a cooler touch and reduced vibration
- Saw Handpiece protected to IP67 for easy cleaning and decontamination
- Power Supply protected to IP65 for use in wet environments
- Includes BD091 Power Supply Box
- Includes BD154 T-Handle Wrench for replacing saw blades

### Compatibility

- Compatible with BD140 Spinal Column Autopsy Saw Blade, 51mm diameter
- Compatible with BD141 Round Autopsy Saw Blade, 64mm diameter
- Compatible with BD142 Section Autopsy Saw Blade, 64mm diameter
- Compatible with BD143 Large Round Autopsy Saw Blade, 76mm diameter



[www.mopec.com](http://www.mopec.com)

800.362.8491

800 Tech Row

Madison Heights, MI 48071