

Hydraulic Cadaver Carrier with False Frame and JZ007 Cadaver Carrier Cover - JA614



SPECIFICATIONS

Construction

- All stainless steel fixed top
- Twin pedestal hydraulics for easy foot-activated elevation
- (4) Heavy-duty 8" casters with total lock
- Bumper guards
- Removable aluminum false frame with JZ007 cadaver cover

Range of Motion

- Elevates from 21.25" to 34" (carrier only)
- Elevates from 37.25" to 50" (with false frame)

Dimensions

- 81" L x 30" W

Controls

- Brake lever on both sides and one on end
- Elevating pedals on both sides

Performance

- Weight capacity: approximately 500 lbs.



www.mopec.com

800.362.8491

800 Tech Row

Madison Heights, MI 48071