

# Mini Formalin Pump - FP025

The Mopec FP025 Mini Formalin Pump is the ultimate portable pumping unit, built for conveniently transferring formalin from Mopec prefilled formalin cubitainers to a dispensing valve. The system is portable and powered by a standard 110v outlet for convenience and ease of use. Customers will find that it meets their needs for quick and convenient dispensing.



## SPECIFICATIONS

### Construction

- All stainless steel construction
- Removable spill containment drip pan with Mopec Formaway Formalin Neutralizing Absorbent Pads

### Dimensions

- 17.1" L x 8.4" W x 15.4" H
- Weight: approximately 23 lbs (10.4 kg)

### Performance

- Demand pump for quick and convenient dispensing
- Portable for convenience and ease of use
- Powered by a standard 110v outlet

### Electrical

- 115v / 1ph / 60hz
- Demand pump: 1.4 GPM, FLA = 0.54A
- 6' IEC320 C13 hospital-grade power cord

### Controls

- Waterproof on/off switch
- Forward sealing dispense valve with flow control

### Compatibility

- Connects to standard cubitainers

### Options & Accessories

- BE092: Mopec Formaway Formalin Neutralizing Absorbent Pads
- BG520: Mopec Prefilled Formalin Cube, 2.5 Gallons
- BG525: Mopec Prefilled Formalin Cube, 5 Gallons



[www.mopec.com](http://www.mopec.com)

800.362.8491

800 Tech Row

Madison Heights, MI 48071