

The Mopec logo features the word "Mopec" in a bold, blue, sans-serif font. A red swoosh underline starts under the 'M' and curves under the 'e'.

Mopec

ELEVATING PATHOLOGY



[MOPEC.COM/AUTOPSY-SAW-BLADES](https://www.mopec.com/AUTOPSY-SAW-BLADES)

Mopec Autopsy Saw Blades

Mopec's affordable autopsy saw blades are hardened for increased longevity and sharpness, earning their title as the highest quality blades available on the market. Our unmatched selection of autopsy saw blades provides you with the options you need for any dissection.

SuperCut Autopsy Saw Blades

Made of the finest tempered stainless steel and coated with an extraordinary "Gold" Titanium Nitride hardening veneer, SuperCut bone cutting blades offer a longer and sharper life than standard autopsy saw blades. Extra sharp teeth provide a clean, precise cut every time.

Compatible with Mopec 810, Mopec 1000 and Stryker 810 Saws

Part ID	Description	Optional Arbor
BD113	Section Blade	With Arbor
BD114	Round Blade	With Arbor

*Please contact a Mopec representative for arborless blades.

Mopec autopsy saw blades are the only blades available with a laser welded arbor to ensure durability and longevity.

Standard Autopsy Saw Blades

Stainless steel autopsy saw blades are hardened to prolong life and sharpness. Mopec's wide selection of stainless steel autopsy blades come with or without arbor and combine the highest quality with an affordable price.

Compatible with Mopec 810, Mopec 1000 and Stryker 810 Saws

Part ID	Description	Optional Arbor
BD101	Large Section Blade	With Arbor
BD104	Round Blade	With Arbor
BD105	Spinal Column Blade	With Arbor
BD108	Small Section Blade	With Arbor

*Please contact a Mopec representative for arborless blades.

Swordfish Autopsy Saw Blades

We've upgraded our Mopec Swordfish saw blades! These new and improved blades are now 10% thicker and made from Stainless Steel instead of Chrome-Plated Carbon Steel!

Compatible with Mopec Swordfish Saw

Part ID	Description	Optional Arbor
BD140	Spinal Column Blade	Without Arbor
BD141	Round Blade	Without Arbor
BD142	Section Blade	Without Arbor
BD143	Large Round Blade	Without Arbor

