

Portable Shelving System

The LD501 portable, five-tier shelving rack conveniently stores within the LD500 cabinet for proper ventilation of its contents. This rack eliminates the need to transfer trays and samples during intra-facility transport. Simply load specimens once for both proper storage and future transport.



SPECIFICATIONS

Construction

- Epoxy-coated posts
- (5) adjustable shelves that can be easily reconfigured (1" increments)
- (2) 5" swivel casters and (2) 5" locking casters made with no metal-exposed parts
- (5) storage trays
- (5) formalin-neutralizing pads for storage trays

Capacity

- Approximate weight per shelf: 200lbs

Dimensions

- Rack: 26.25" L x 26.25" W x 67.5" H
- Storage Trays: 22.5" L x 17.5" W x 3" H

Performance

- During transport, smaller containers are conveniently located in trays for spill-proof containment while absorbing and neutralizing formalin

Compatibility

- LD500 Ventilated Tissue Cabinet

Options & Accessories

- JZ011: Optional transport cover available

Related/Alternate/Obsolete Product Model Numbers

Model No.	Description	Status	Notes
SC614-S	Portable Shelving System	Replaced	Use LD501

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI