


Pedestal Autopsy Table, Vented, with Optional Elevation

1036-41-ATD



1036-41-TD



SPECIFICATIONS

Construction

- Heliarc-welded seams and joints, ground and polished to match adjoining surfaces
- (4) Fully perforated grid plates
- Mixing faucet with gooseneck spout and wrist blade handles
- Heavy duty waste disposal (1/2 HP)
- Deluxe hydro-aspirator with reverse flow, cold water control valve with 10ft of hose
- Spray hose assembly with cold water control valve, nozzle with 10ft of hose

Materials

- 304 stainless steel construction with a polished No. 4 finish

Utilities / Power

- (1) 115v / 60hz / 1ph / 15a circuit for station
- (1) 115v / 60hz / 1ph / 15a circuit for disposer
- Exhaust recommendation: 450-500 CFM

Dimensions

- 1036-41-TD: 104"L x 32"W x 35"H
- 1036-41-ATD: 104"L x 32"W x 32" - 44"H
- Sink depth: 8"

Capacity

- 750 lbs

Performance

- GFCI protected receptacles with waterproof covers
- Back-flow protection for all plumbing fixtures to prevent back-flow of contaminated water
- Rapid and positive drainage
- (2) Access panels on pedestal

Options & Accessories

- Elevation Optional: 1036-41-ATD
- CO017: Scale stand, table mounted
- CO009: Perimeter rinse
- CO012: In/Cm Scale/Ruler

Alternate Product Model Numbers

Model No.	Description	Status	Notes
CE100	Pedestal Autopsy Table, Vented	Replaced	Use 1036-41-TD
CE200	Pedestal Autopsy Table, Vented, Elevating	Replaced	Use 1036-41-ATD

Related / Obsolete Product Model Numbers

- CE110, CE210, CE340, CE350, CE400, CE500, CE800, CE850, CE900
- 1036-41, 1036-41-A, 1036-41-SR, 1036-41-MR, 1036-24, 1036-28

Variations

Model No.	Description	Dimensions
1036-41-TD	Pedestal Autopsy Table, Vented	104" L x 32" W x 35" H
1036-41-ATD	Pedestal Autopsy Table, Vented, Elevating	104" L x 32" W x 32" - 44" H

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI