


Dissection Table, Dip Tank with Hand Crank Elevation

This Mopec dissection table offers unmatched versatility and functionality for anatomy and autopsy professionals. Its integrated dip tank ensures efficient specimen collection and preservation, while the manual hand crank elevation allows precise adjustments for ergonomic comfort during procedures.



SPECIFICATIONS

Construction

- Ergonomically placed hand crank for easier and safer elevation of internal tray
- Magnetic latches to secure cover lids
- Holes in tray for a leveled specimen during elevation
- Flanged drain hole in tray
- (4) heavy duty 5" casters for easier transporting

Materials

- 304 stainless steel construction with a polished No. 4 finish

Dimensions

- 96" L x 32.6" W x 38" H
- Interior width for cadaver: 27"

Capacity

- 300 Lbs
- Max fluid capacity: 119.7 gallons

Options & Accessories

- HO005: Post leg-holders
- HO006: Book holder
- HO030: Instrument tray holder

Related/Alternate/Obsolete Product Model Numbers

Model No.	Description	Status	Notes
600022-01	6" Single Connection	Obsolete	
600022-02	6" Single Connection	Obsolete	
600022-03	6" Single Connection	Obsolete	
600022-04	6" Single Connection	Obsolete	

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI