

Bridge-Style Leg Rest and Holder for Dissection Table



SPECIFICATIONS

Construction

- Built to support both legs of a cadaver at the same height and width

Notes

- Must remove one grid plate to use on hooded HB series carts
- Hoods will not close when using on HB series carts

Materials

- 304 stainless steel construction with a polished No. 4 finish

Dimensions

- 26.7" L x 10" W x 17" H

Compatibility

- HA100
- HA400
- HB1400
- HB1410
- HB1500
- HB1510
- HB1550

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI