

Cadaver Cover, Vinyl, Heavy Duty

Allows for the discrete transport of a patient from the hospital room to the mortuary. The rectangular shape of the cover conceals the body giving the illusion of a standard gurney while hiding the cadaver. Cadaver covers are made of black out material which can be easily cleaned.



SPECIFICATIONS

Materials

- Heavy duty vinyl type material
- Color: white or blue

Dimensions

- Vary depending on option

Options:

- JZ006: 81" L x 31" W x 18" D, white
- JZ007: 81" L x 31" W x 18" D, blue
- JZ015: 81" L x 33" W x 21" D, blue
- JZ016: 81" L x 33" W x 21" D, white
- JZ020: 81" L x 27" W x 18" D, white
- JZ021: 81" L x 27" W x 18" D, blue

CLEANING INSTRUCTIONS

1. This material is not machine washable.
2. Use a very mild soap and water solution and wipe with a damp cloth.
3. Air dry or use a separate dry cloth to wipe the material dry.

Features

Model No.	Material	Color	Dimensions	Compatibility
JZ006	Heavy Duty Vinyl	White	81" L x 31" W x 18"H	JA600
JZ007	Heavy Duty Vinyl	Blue	81" L x 31" W x 18" H	JA600
JZ015	Heavy Duty Vinyl	Blue	81" L x 33" W x 21"H	GA110
JZ016	Heavy Duty Vinyl	White	81" L x 33" W x 21" H	GA110
JZ020	Heavy Duty Vinyl	White	81" L x 27" W x 18"H	JA611 / JA630
JZ021	Heavy Duty Vinyl	Blue	81" L x 27" W x 18" H	JA611 / JA630



www.mopec.com

800.362.8491

800 Tech Row

Madison Heights, MI 48071