


Ensemble Pegboard Requisition Form Holder



SPECIFICATIONS

Construction

- Designed for the Mopec Maestro's Ensemble Pegboard
- Designed to fit 8 1/2" x 11" sized requisitions with room to spare

Materials

- 304 stainless steel construction with a polished No. 4 finish

Dimensions

- 10.3" W x 9.0" H x 1.1" D
- Internal: 10.25" W x 8.25" H x 1.0" D

Compatibility

- Fits all Mopec Pegboards, including AX031R and AX031L pegboard wings

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI