

Ensemble Pegboard Shelves



AX006



AX072

SPECIFICATIONS

Construction

- Designed for the Mopec Maestro's Ensemble Pegboard

Materials

- 304 stainless steel construction with a polished No. 4 finish

Dimensions

- AX006: 14.0" L x 8.0" D x 2.0" H
- AX072: 7.0" L x 8.0" D x 2.0" H

Capacity

- AX006: 12 lbs
- AX072: 12 lbs

Compatibility

- Fits all Mopec Pegboards, including AX031R and AX031L pegboard wings
- Can also be attached to the front of grossing stations, above the vent grille

Variations

Model No.	Description	Dimensions	Capacity
AX006	Pegboard Shelf, 14"	14" L x 8" W x 2" H	12 lbs
AX072	Pegboard Shelf, 7"	7" L x 8" W x 2" H	12 lbs

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI