


Mopec Maestro PathCam Kit for 2018-2022 Maestro Platforms

Upgrade your existing grossing station with a PathCam kit to transform your workflow and bring your lab into the digital age. These kits provide seamless integration, enabling advanced imaging and annotation capabilities without replacing your entire workstation. Designed for 2018–2022 Maestro platforms.



SPECIFICATIONS

Construction

- (1) 15a receptacle, 125 vac, 35mm din rail mount to supply camera power
- (1) Panel-Mount USB-A 2.0 Coupler with cover plate
- Includes all necessary hardware to complete installation

Materials

- Bracket: 304 stainless steel construction

Dimensions

- Mounting plate: 5.8" L x 8.8" W x 4.6" H
- Receptacle: 5.7" L x 2.1" W x 2.3" H
- Panel-Mount USB with cover plate: 1.7" L x 2.0" W x 4.3" H

Compatibility

- 2018–2022 Maestro platform
- AX265HD camera
- AX266HR camera

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI