


NOTE:
PARTS ARE FOR
KITING PURPOSES ONLY


CAMERA IS REFERENCE ONLY

REV	DESCRIPTION	DATE	ENR	APP
A	INITIAL RELEASE	28-JUN-22		JJ
				DA
				JJ

DETAIL A



NOTICE LOCATION CUT MARKS FOR HOLSTER



Item Number	Quantity	File Name	Title	Material	Rev
1	1	M-4311.psm	USB PANEL MOUNT BRACKET	SHT, SS, 18 GA, 304	A
2	2	PB0016.par	NUT - SS, 1/4"-20, NYLON LOCK	McMASTER-CARR # 91831A029	A
3	2	PB0247.par	1/4"-20 x 1/2" STUD SS		A
4	4	PB0338.par	#6-32 INSERT SPIN TITE	MARSHALL SALES # AE-KS-832-80	A
5	1	PB0512.par	SCREW - SS, 10-32 x 3/8", BUTTON HEAD	McMASTER-CARR # 92949A263	A
6	1	PE0134.par	LINE CORD 12/3		
7	1	PE0954.par	POWER KIT, PATHCAM FOR STANDARD ENCORE		A
8	1	PE0955.par	SD CARD, PATHCAM ENCORE		A
9	1	PE0957.par			
10	1	S6883-F.psm	CAMERA HOLSTER	SHT, SS 11 GA, 304	E
11	1	S7586-B.psm	MAESTRO, PATHCAM COVER PLATE, FOR AXI23	SHT, SS 11 GA, 304	B

THIS DOCUMENT CONTAINS PROPRIETARY INFORMATION OF MOPEC, AND MAY NOT BE COPIED IN WHOLE OR PART, AND NO USE SHALL BE MADE OF IT OTHER THAN FOR WHICH IT HAS BEEN SUBMITTED WITHOUT THE WRITTEN PERMISSION OF A MOPEC MANAGER	UNLESS OTHERWISE SPECIFIED ALL DIMENSIONS IN INCHES TOLERANCES 0.00 ± 0.1 in 0.000 ± 0.05 in 0.0000 ± 0.001 in 0.00000 ± 0.0010 in ANGLES: ±0.5° BREAK ALL CORNERS 0.063 x 45° ALL FILLET RADI R0.063 (0.125in) INSPECT OVAL DIMENSIONS	MATERIAL: <
--	---	--