


Mopec Maestro PathCam Kit for 2023-2024 Maestro Platforms

Upgrade your existing grossing station with a PathCam kit to transform your workflow and bring your lab into the digital age. These kits provide seamless integration, enabling advanced imaging and annotation capabilities without replacing your entire workstation. Designed for 2023-2024 Maestro platforms.



*Camera not included



SPECIFICATIONS

Construction

- (1) Panel-Mount USB-A 2.0 Coupler with cover plate
- Includes all necessary hardware to complete installation

Materials

- Bracket: 304 stainless steel construction

Dimensions

- Mounting plate: 5.8" L x 8.8" W x 4.6" H
- Receptacle: 5.7" L x 2.1" W x 2.3" H
- Panel-Mount USB with cover plate: 1.7" L x 2.0" W x 4.3" H

Compatibility

- 2023-2024 Maestro platform
- AX265HD camera
- AX266HR camera

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI