

NOTE:
PARTS ARE FOR
KITING PURPOSES ONLY

CAMERA IS REFERENCE ONLY

REV	DESCRIPTION	DATE	ENR	APP
A	INITIAL RELEASE	28-JUN-22		JJ

DETAIL A

ATTACH ONTO DIN RAIL,
WHERE POWER SUPPLY IS LOCATED
WIRE TO TO 110V DIN RAIL
CONNECTORS IN SIDE PANEL

NOTICE LOCATION CUT MARKS FOR HOLSTER

Item #	Qty	File Name	Title	Material	Rev
		AX043-L3.asm	MAESTRO, PATHCAM SYSTEM, ADVANCED		A
1	1	M-431l.psm	USB PANEL MOUNT BRACKET	SHT, SS, 1/8 GA, 304	A
2	2	PB0016.par	NUT - SS, 1/4"-20, NYLON LOCK	McMASTER-CARR # 91831A029	A
3	2	PB0247.par	1/4"-20 x 1/2" STUD SS		A
4	4	PB0338.par	#8-32 INSERT SPIN TITE	MARSHALL SALES # AE-KS-832-80	A
5	1	PB0512.par	SCREW - SS, 10-32 x 3/8", BUTTON HEAD	McMASTER-CARR # 92949A263	A
6	1	PB0780.par	2-3/4" Carbide Tipped Hole Cutter, 1" Depth of Cut		A
7	1	PB0785.asm			
7.1	1	4289A35-Changeable-Bi-ade Deburring Scraper.par			
8	1	PE0134.par	LINE CORD 12/3		
9	1	PE0754.par	RECEPTACLE, 15A 125 VAC, 35MM DIN RAIL MOUNT	AUTOMATION DIRECT FA-REC2	A
10	1	PE0954.par	POWER KIT, PATHCAM FOR STANDARD ENCORE		A
11	1	PE0955.par	SD CARD, PATHCAM ENCORE		A
12	1	PE0957.asm			
13	1	PM0733.par	3M DOUBLE SIDED TAPE 1" WIDE		A
14	1	S6883-F.psm	CAMERA HOLSTER	SHT, SS 11 GA, 304	E
15	1	S7586-B.psm	MAESTRO, PATHCAM COVER PLATE, FOR AXI23	SHT, SS 11 GA, 304	B

THIS DOCUMENT CONTAINS
PROPRIETARY INFORMATION OF
MOPEC, AND MAY NOT BE COPIED IN WHOLE
OR PART, AND NO USE SHALL BE MADE OF
IT OTHER THAN FOR WHICH IT HAS BEEN
SUBMITTED WITHOUT THE WRITTEN
PERMISSION OF A MOPEC MANAGER

MOPEC 2021

UNLESS OTHERWISE SPECIFIED
ALL DIMENSIONS IN INCHES
TOLERANCES
0.0 ± 0.1 in
0.00 ± 0.05 in
0.000 ± 0.001 in
0.0000 ± 0.0010 in
ANGLES
± 0.5°
BREAK ALL CORNERS 0.063 x 45°
ALL FILLET RADI R0.063 (0.025 in)
RESPECT OVAL DIMENSIONS

MATERIAL:



DRAWN BY: J. Shanks
6/10/2024
CHECKED BY:
DATE: 6/10/2024
APPROVED BY:
DATE: MMCLAIN 6/10/2024

PART NUMBER	AX043-L3
DESCRIPTION	MAESTRO, PATHCAM SYSTEM, ADVANCED
SIZE	A2
FILE NAME	AX043-L3.dft
SCALE	1:1.25
DO NOT SCALE IF IN DOUBT ASK	
SHEET 1 OF 2	

AX043-L3.dft