

Mopec PathCam Kit for Legacy MB, ML, ME, MG Stations

Upgrade your existing grossing station with a PathCam kit to transform your workflow and bring your lab into the digital age. These kits provide seamless integration, enabling advanced imaging and annotation capabilities without replacing your entire workstation. Designed for MB, ML, ME, and MG platforms.



*Camera not included



SPECIFICATIONS

Construction

- (1) 15a receptacle, 125 vac, 35mm din rail mount to supply camera power
- (1) Panel-Mount USB-A 2.0 Coupler with cover plate
- (1) 2-3/4" hole saw included to core hole for camera lens
- (1) Deburring Tool for camera cutout
- Includes all necessary hardware to complete installation

Materials

- Bracket: 304 stainless steel construction

Dimensions

- Mounting plate: 5.8" L x 8.8" W x 4.6" H
- Receptacle: 5.7" L x 2.1" W x 2.3" H
- Panel-Mount USB with cover plate: 1.7" L x 2.0" W x 4.3" H

Compatibility

- Legacy Stations: MB, ML, ME, MG
- AX265HD camera
- AX266HR camera

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI