


AX072-B.dft



REV	DESCRIPTION	DATE	ECN #	APPR
A	INITIAL RELEASE	00-XXX-19	XXXXX	MM
B	INVERTED	3/29/21	1038	JJ

Item Number	Quantity	File Name	Title	Material	Revision
1	1	S6800-Apsm	PEGBOARD MEDIUM MOUNTING TAB	18 GA. TYPE T304 SS	A
2	1	S6880.psm	ADJUSTABLE SHELF/PEGBOARD	18 GA. TYPE T304 SS	B

THIS DOCUMENT CONTAINS PROPRIETARY INFORMATION OF MOPEC, AND MAY NOT BE COPIED IN WHOLE OR PART, AND NO USE SHALL BE MADE OF IT OTHER THAN FOR WHICH IT HAS BEEN SUBMITTED WITHOUT THE WRITTEN PERMISSION OF A MOPEC MANAGER

MOPEC 2024

UNLESS OTHERWISE SPECIFIED:
ALL DIMENSIONS IN INCHES

TOLERANCES

0.0	± 0.1 in
0.00	± 0.05 in
0.000	± 0.010 in
0.0000	± 0.0010 in
ANGLES	± 0.5°

BREAK ALL CORNERS 0.063 X 45°
ALL FILLET RADII R0.063 /0.125in
INSPECT OVAL DIMENSIONS

MOPEC GROUP

MOPEC MORTTECH

ILLUSTRATION TECHNOLOGY MANUFACTURING

DRAWN TO 3RD ANGLE PROJECTION

3RD ANGLE PROJECTION

DRAWN BY: B. LAY
DATE: 12/5/2024

CHECKED BY: ---
DATE: 12/5/2024

APPROVED BY: M. McCLAIN
DATE: 12/5/2024

PART NUMBER: AX072-B

DESCRIPTION: MAESTRO, PEGBOARD, 7IN SHELF

SIZE: A2

FILE NAME: AX072-B.dft

SCALE: 1:1

DO NOT SCALE IF IN DOUBT ASK

SHEET 2 OF 2