

## DISSECTION TABLES

# Covered Dissection Table with Lid Assist, Elevating, Optional Ventilation

The Mopec Covered Dissection Tables combine safety, functionality, and discretion to meet the needs of anatomy labs and teaching facilities. Featuring the exclusive and secure Mopec lid assist, the hinged covers lock neatly under the table during use, minimizing disruptions and enhancing safety. Elevation ensures ergonomic comfort, making this table an ideal choice for fostering a safe and efficient dissection experience.



## SPECIFICATIONS

### Construction

- All stainless-steel construction
- Heliarc-welded seams and joints, ground and polished to match adjoining surfaces
- Built-in tapered slope
- Includes drain hole with PVC drain valve
- (4) heavy duty 8" casters with (3) total lock and (1) directional lock
- (4) Post sockets included

### Materials

- 304 stainless steel construction with a polished No. 4 finish

### Dimensions

- 92" L x 30" W x 32" – 39.9" H

### Capacity

- 350 Lbs

### Performance

- Mopec Lid Assist: hinged cover with swing latches for easy open and close
- Hand crank elevation: 32" – 39.9"
- Exhaust requirements: 220–340 CFM (HB1510)
- Static pressure: .34–.67" WG (HB1510)
- ADA wheelchair accessible (HB1550)

### Options & Accessories

- Ventilation: Use HB1510
- Non-Elevation: See HB1400 series
- AX107: Bucket holder
- HO001: Bridge-style leg holder
- HO002: Bridge-style book holder
- HO005: Post leg-holders
- HO006: Post book holder
- HO024: Full-length shelf
- HO030: Instrument tray holder
- HO032: Short shelf (use with AX107)
- HO1110: Misting system

## Related/Alternate/Obsolete Product Model Numbers

Model No.	Description	Status	Notes
600021-HC-LA	Dissection Cart, Elevated, Lid Assist	Replaced	Use HB1400
600024-HC-LA	Dissection Cart, Ventilated, Elevated, Lid Assist	Replaced	Use HB1410
HB1520	Dissection Cart, Dual Ventilated, Elevated, Lid Assist	Obsolete	

## Variations

Model No.	Description	Ventilation	Elevation
HB1500	Dissection Cart with Lid Assist	Non-Ventilated	Elevating
HB1510	Dissection Cart with Lid Assist	Single Duct Ventilation	Elevating
HB1550	Dissection Cart with Lid Assist	Non-Ventilated	Elevating, ADA

**MOPEC**  
GROUP

**mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI