


Cadaver Cover, Vinyl, Standard Duty

Allows for the discrete transport of a patient from the hospital room to the mortuary. The rectangular shape of the cover conceals the body giving the illusion of a standard gurney while hiding the cadaver. Cadaver covers are made of black out material.



SPECIFICATIONS

Materials

- Rubber-backed blackout fabric
- Color: white

Dimensions

- Vary depending on option

Options:

- JZ005: 81" L x 31" W x 18" D
- JZ010: 81" L x 40" W x 20" D
- JZ014: 81" L x 33" W x 21" D
- JZ017: 81" L x 27" W x 18" D
- JZ018: 81" L x 31.62" W x 18" D
- JZ019: 87" L x 31" W x 18" D

CLEANING INSTRUCTIONS

1. Choose a very mild detergent and follow the instructions. Do not use bleach.
2. Wash in cold water on the gentle cycle.
3. Air dry or dry on low heat.
4. Wash as needed.

Features

Model No.	Material	Color	Dimensions	Compatibility
JZ005	Vinyl	White	81" L x 31" W x 18"H	JA600
JZ010	Vinyl	White	81" L x 40" W x 20" H	DC410
JZ014	Vinyl	White	81" L x 33" W x 21"H	GA110
JZ017	Vinyl	White	81" L x 27" W x 18" H	JA611 / JA630
JZ018	Vinyl	White	81" L x 31.62" W x 18"H	JA612 / JA631
JZ019	Vinyl	White	87" L x 31" W x 18" H	HB300



www.mopec.com

800.362.8491

800 Tech Row

Madison Heights, MI 48071