


Ensemble Pegboard Glasses Holster



SPECIFICATIONS

Construction

- Designed for the Mopec Maestro's Ensemble Pegboard
- Designed to hold your PPE and safety glasses

Capacity

- 2.5 lbs

Materials

- 304 stainless steel construction with a polished No. 4 finish

Dimensions

- 3.8" L x 5.6" D x 0.8" H

Compatibility

- Fits all Mopec Pegboards, including AX031R and AX031L pegboard wings
- Can also be attached to the front of grossing stations, above the vent grille

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI