


Ensemble Hitch Writing Ledge



SPECIFICATIONS

Construction

- Mounts onto “trailer hitch” that attaches onto the front of the station, below the work surface

Capacity

- 5 lbs

Materials

- 304 stainless steel construction with a polished No. 4 finish

Performance

- Can go on either the left or right side of the station
- Can be placed onto “trailer hitch” mount when in use and either pivoted out of the way or removed when not in use

Dimensions

- 11.5” L x 11.6” D x 3.3” H

Compatibility

- Fits all Mopec Maestro and Maestro Encore grossing stations

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI