


# Ensemble Light Bar



## SPECIFICATIONS

### Construction

- Designed to mount under the shelf/above the vent grille on the Maestro Encore Grossing station
- Daylight white light, 6500K

### Materials

- Painted Metal frame with a polycarbonate lens, LED light diodes

### Dimensions

- 24" L x 2.2" D x 1.3" H

### Compatibility

- Fits all Mopec Maestro and Maestro Encore grossing stations

## NOTE:

- This option **must** be ordered and at the time of purchase and installed during manufacturing of the grossing station. It cannot be installed in the field after delivery.

**MOPEC**  
GROUP

**Mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI