

Exploded View Drawing of a Dissection Board Assembly

Grid Labels: A, B, C, D, E, F, G, H (Vertical); 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 11 (Horizontal)

Callouts: 1/8, 2/2

Parts List:

Item Number	Quantity	File Name	Title	Material	Revision
		AX064-096-A.asm	MAESTRO, DISSECTION BOARD, WHITE, 54"		A
1	8	PB0609.par	TAPPING INSERT, 1/4"-20 INTERNAL, BRASS	Stainless steel	A
2	2	S7540-A.par	DISSECTION BOARD, 21" x 18.5", WHITE	3/4" HDPE - NATURAL, TEXTURED	A

Notes:

- THIS DOCUMENT CONTAINS PROPRIETARY INFORMATION OF MOPEC, AND MAY NOT BE COPIED IN WHOLE OR PART, AND NO USE SHALL BE MADE OF IT OTHER THAN FOR WHICH IT HAS BEEN SUBMITTED WITHOUT THE WRITTEN PERMISSION OF A MOPEC MANAGER
- UNLESS OTHERWISE SPECIFIED: ALL DIMENSIONS IN INCHES
TOLERANCES: 0.0 ± 0.1 in, 0.00 ± 0.05 in, 0.000 ± 0.010 in, 0.0000 ± 0.0010 in, ANGLES ±0.5°
- BREAK ALL CORNERS 0.063 X 45°
ALL FILLET RADII R0.063 /0.125in
INSPECT OVAL DIMENSIONS


Company Information: MOPEC GROUP, MOPEC MORTTECH, MANUFACTURING

Drawn To: 3RD ANGLE PROJECTION


Drawn By: A. SALAME
Date: 1/28/2025
Checked By: ---
Date: 1/28/2025
Approved By: MMCLAIN
Date: 1/28/2025

Part Number: AX064-096-A
Description: MAESTRO, DISSECTION BOARD, WHITE, 54"
Size: A2
File Name: AX064-096-A.dft
Rev: A
Scale: 1:2
Do Not Scale If In Doubt Ask
Sheet: 1 OF 3

AX064-096-A.dft



AX064-096-A.dft



CUTTING BOARD IS TO
BE FREE FROM SCRATCHES,
STAINS, DISCOLORATION, &
DEFORMATION

REV	DESCRIPTION	DATE	ECN #	APPR
A	INITIAL RELEASE	00-XXX-19	XXXXX	MM

TOP

BOTTOM


THIS DOCUMENT CONTAINS
PROPRIETARY INFORMATION OF
MOPEC, AND MAY NOT BE COPIED IN WHOLE
OR PART, AND NO USE SHALL BE MADE OF
IT OTHER THAN FOR WHICH IT HAS BEEN
SUBMITTED WITHOUT THE WRITTEN
PERMISSION OF A MOPEC. MANAGER

MOPEC 2024

UNLESS OTHERWISE SPECIFIED:
ALL DIMENSIONS IN INCHES

TOLERANCES	0.0	± 0.1 in
	0.00	± 0.05 in
	0.000	± 0.010 in
	0.0000	± 0.0010 in
	ANGLES	±0.5°

BREAK ALL CORNERS 0.063 X 45°
ALL FILLET RADII R0.063 /0.125in
INSPECT OVAL DIMENSIONS



MOPEC
GROUP



MORTECH
MANUFACTURING

DRAWN TO
3RD ANGLE
PROJECTION



DRAWN BY:
A. SALAME
DATE: 1/28/2025

CHECKED BY:
DATE: ---
1/28/2025

APPROVED BY:
DATE: M.McCLAIN
1/28/2025

PART
NUMBER: AX064-096-A

DESCRIPTION: MAESTRO, DISSECTION BOARD, WHITE, 54"

SIZE: A2	FILE NAME: AX064-096-A.dft	REV: A
SCALE: 1:2.67	DO NOT SCALE IF IN DOUBT ASK	SHEET 3 OF 3