

Ensemble Shelf-Mounted Dual Gooseneck LED Task Lights



SPECIFICATIONS

Construction

- Mount to the shelf area of the grossing station, evenly spaced apart, one for each half of the work surface

Materials

- ABS plastic, rubber coated flexible metal arm, stainless steel hardware

Dimensions

- 6.5" L x 2.4" W x 17.3" H each
- Flexible arm: 15.3"
- Flexible arm with light head: 16.5"

Compatibility

- Fits all Mopec Maestro and Maestro Encore grossing stations

Options & Accessories

- For one light , please order AX081 Single Gooseneck LED Task Light

NOTE:

- We recommend ordering at the time of purchase to be properly installed during manufacturing of the grossing station.

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI