

PathCam Monitor, 21" Touchscreen



SPECIFICATIONS

Construction

- Used to display PathCam Imaging Systems

Utilities / Power

- 115 V

Dimensions

- 21" display screen

Performance

- LED full HD display
- 1920 x 1080 pixels
- 10-point multi-touch interactivity
- Inputs: DisplayPort, VGA, and HDMI
- (2) integrated speakers
- VESA mount compatible
- Blue light reduction
- HID compliant

Compatibility

- Maestro and Maestro Encore grossing stations equipped with PathCam Imaging Systems

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI