


# PathCam Hands-Free Foot Pedal Command



## SPECIFICATIONS

### Construction

- Designed to allow the user to program four different pedals to the functions of their choosing
- Includes 2 meter USB cable

### Utilities / Power

- Self-powered from USB port

### Dimensions

- 16.6" L x 5.25" W x 1.25" H

### Controls

- Programmable quad USB foot switch
- Side-by-side footswitch configuration

### Performance

- Simple plug-and-play USB compliant foot switch which uses standard USB HID driver

### Compatibility

- AX265 PathCamHD Camera
- AX266 PathCamHR Camera
- AX102 PathCam Micron Microscope

**MOPEC**  
GROUP

**mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI