

Maestro Corner Dolly, Set of 4



SPECIFICATIONS

Construction

- Designed to sit under the 4 corners of a Maestro or Maestro Encore grossing station to assist Mopec installation personnel while transporting

Materials

- 304 stainless steel
- Non-marking rubber caster wheels

Dimensions

- 14.1" L x 14.1" W x 4.8" H

Capacity

- 800 lbs

Compatibility

- Fits the Mopec Maestro and Maestro Encore grossing stations

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI