

REV	DESCRIPTION	DATE	ECN #	APPR
A	INITIAL RELEASE	00-XXX-19	XXXXX	MM

THIS DOCUMENT CONTAINS PROPRIETARY INFORMATION OF MOPEC, AND MAY NOT BE COPIED IN WHOLE OR PART, AND NO USE SHALL BE MADE OF IT OTHER THAN FOR WHICH IT HAS BEEN SUBMITTED WITHOUT THE WRITTEN PERMISSION OF A MOPEC. MANAGER

MOPEC 2024

UNLESS OTHERWISE SPECIFIED:
ALL DIMENSIONS IN INCHES

TOLERANCES

0.0	± 0.1 in
0.00	± 0.05 in
0.000	± 0.010 in
0.0000	± 0.0010 in
ANGLES	± 0.5°

BREAK ALL CORNERS 0.063 X 45°
ALL FILLET RADII R0.063 /0.125in
INSPECT OVAL DIMENSIONS

MOPEC GROUP

MORTECH

DRAWN TO 3RD ANGLE PROJECTION

DRAWN BY: B. LAY
DATE: 8/22/2022

CHECKED BY: ---
DATE: 8/22/2022

APPROVED BY: M. McClain
DATE: 8/22/2022

PART NUMBER: AX117L-A

DESCRIPTION: ARM MOUNTING BRKT, SCREEN SIDE

SIZE: A2

SCALE: 1:5


DO NOT SCALE IF IN DOUBT ASK

SHEET 1 OF 2

REV: A

AX117L-A.dft

REV	DESCRIPTION	DATE	ECN #	APPR
A	INITIAL RELEASE	00-XXX-19	XXXXXX	MM



Item Number	Quantity	File Name	Title	Material	Rev
1	1	M-4121-Apsm	MB, HOOD, PILLAR SLIDE	18 GA. TYPE T304 SS	A
2	4	S2046-Apsm	GUSSET	11 GA. TYPE T304 SS	A
3	1	S7374-Apsm	MAESTRO, TRACK ,MOUNTING BRKT	14 GA. TYPE T304 SS	A
4	1	S7487-Apsm	END-CAP MOUNTING BRKT	18 GA. TYPE T304 SS	A

THIS DOCUMENT CONTAINS
PROPRIETARY INFORMATION OF
MOPEC, AND MAY NOT BE COPIED IN WHOLE
OR PART, AND NO USE SHALL BE MADE OF
IT OTHER THAN FOR WHICH IT HAS BEEN
SUBMITTED WITHOUT THE WRITTEN
PERMISSION OF A MOPEC MANAGER

UNLESS OTHERWISE SPECIFIED:
ALL DIMENSIONS IN INCHES
TOLERANCES

0.0	± 0.1 in
0.00	± 0.05 in
0.000	± 0.010 in
0.0000	± 0.0010 in
ANGLES	±0.5°

BREAK ALL CORNERS 0.063 X 45°
ALL FILLET RADII R0.063 /0.125in
INSPECT OVAL DIMENSIONS



MOPEC
GROUP

mopec **MORTECH**
ELEVATING PATHOLOGY MANUFACTURING

DRAWN TO
3RD ANGLE
PROJECTION

DRAWN BY:
B. LAY
DATE: 8/22/2022

CHECKED BY:

DATE: 8/22/2022

APPROVED BY:
M.McCLAIN
DATE: 8/22/2022

PART NUMBER		AXI17L-A	
DESCRIPTION: ARM MOUNTING BRKT, SCREEN SIDE			
SIZE: A2	FILE NAME: AXI17L-A.dft		REV: A
SCALE: 1:2	DO NOT SCALE IF IN DOUBT ASK	SHEET 2 OF 2	

AX17L-A.dft