




REV	DESCRIPTION	DATE	ECN #	APPR
A	INITIAL RELEASE	00-XXX-19	XXXXX	MM

Item #	Qty.	File Name	Title	Material	Rev.
1'	1	M-4121-A	MB, HOOD, PILLAR SLIDE	18 GA. TYPE T304 SS	A
2'	6	PB0039	SCREW - SS, 8-32 x 1/2"L, TRUSS PHILLIPS	McMASTER-CARR # 91770A194	
3	6	PB0338	#8-32 INSERT SPIN TITE	Stainless steel	A
4	4	S2046-A	GUSSET	11 GA. TYPE T304 SS	A
5	1	S7374-A	MAESTRO, TRACK MOUNTING BRKT	14 GA. TYPE T304 SS	A
6'	1	S7487-A	END-CAP MOUNTING BRKT	18 GA. TYPE T304 SS	A

THIS DOCUMENT CONTAINS PROPRIETARY INFORMATION OF MOPEC, AND MAY NOT BE COPIED IN WHOLE OR PART, AND NO USE SHALL BE MADE OF IT OTHER THAN FOR WHICH IT HAS BEEN SUBMITTED WITHOUT THE WRITTEN PERMISSION OF A MOPEC MANAGER

MOPEC 2024

UNLESS OTHERWISE SPECIFIED:
ALL DIMENSIONS IN INCHES

TOLERANCES	0.0	± 0.1 in
	0.00	± 0.05 in
	0.000	± 0.010 in
	0.0000	± 0.0010 in
ANGLES		±0.5°

BREAK ALL CORNERS 0.063 X 45°
ALL FILLET RADII R0.063 /0.125in
INSPECT OVAL DIMENSIONS

MOPEC GROUP

MOPEC MORTTECH

subcontractor MANUFACTURING

DRAWN TO 3RD ANGLE PROJECTION

DRAWN BY: B. LAY
DATE: 9/14/2022

CHECKED BY: ---
DATE: 9/14/2022

APPROVED BY: M. McCLAIN
DATE: 9/14/2022

PART NUMBER: AX117R-A

DESCRIPTION: ARM MOUNTING BRKT, SCREEN SIDE

SIZE: A2

FILE NAME: AX117R-A.dft


SCALE: 1:2

DO NOT SCALE IF IN DOUBT ASK

SHEET 1 OF 2

REV: A

AX117R-A.dft



REV	DESCRIPTION	DATE	ECN #	APPR
A	INITIAL RELEASE	00-XXX-19	XXXXX	MM

Item Number	Quantity	File Name	Title	Material	Rev
1"	1	M-4121-Apsm	MB, HOOD, PILLAR SLIDE	18 GA. TYPE T304 SS	A
2"	6	PB0039.par	SCREW - SS, 8-32 x 1/2"L, TRUSS PHILLIPS	McMASTER-CARR # 91770A194	
3	6	PB0338.par	#8-32 INSERT SPIN TITE	Stainless steel	A
4	4	S2046-Apsm	GUSSET	11 GA. TYPE T304 SS	A
5	1	S7374-Apsm	MAESTRO, TRACK, MOUNTING BRKT	14 GA. TYPE T304 SS	A
6"	1	S7487-Apsm	END-CAP MOUNTING BRKT	18 GA. TYPE T304 SS	A

THIS DOCUMENT CONTAINS PROPRIETARY INFORMATION OF MOPEC, AND MAY NOT BE COPIED IN WHOLE OR PART, AND NO USE SHALL BE MADE OF IT OTHER THAN FOR WHICH IT HAS BEEN SUBMITTED WITHOUT THE WRITTEN PERMISSION OF A MOPEC. MANAGER

MOPEC 2024

UNLESS OTHERWISE SPECIFIED:
ALL DIMENSIONS IN INCHES

TOLERANCES

0.0	± 0.1 in
0.00	± 0.05 in
0.000	± 0.010 in
0.0000	± 0.0010 in
ANGLES	±0.5°

BREAK ALL CORNERS 0.063 X 45°
ALL FILLET RADII R0.063 /0.125in
INSPECT OVAL DIMENSIONS

DRAWN TO 3RD ANGLE PROJECTION

MOPEC GROUP

MOPEC MORTECH
Subcontractor MANUFACTURING

DRAWN BY: B. LAY
DATE: 9/14/2022

CHECKED BY: ---
DATE: 9/14/2022

APPROVED BY: M.McCLAIN
DATE: 9/14/2022

PART NUMBER: AXI17R-A

DESCRIPTION: ARM MOUNTING BRKT, SCREEN SIDE

SIZE: **A2** FILE NAME: AXI17R-A.dft REV: **A**

SCALE: **1:1** DO NOT SCALE IF IN DOUBT ASK

SHEET 2 OF 2

AX117R-A.dft