

# Maestro Encore Safety Sash

Offering critical protection from hazardous formalin fumes, chemical spills and splashes, the Safety Sash ensures a safer laboratory environment without obstructing visibility or workflow. With seamless integration for Ensemble pegboard accessories, formalin dispensing, and the PathCam Imaging System, the AX124 Safety Sash is an essential tool for maximizing safety and efficiency in any pathology lab.



AX124L-048



AX124R-048



AX124L-060



AX124R-060



AX124L-072



AX124R-072

## SPECIFICATIONS

### Construction

- Designed to reduce exposure to formalin fumes and hazardous chemical spills and splashes
- Engineered to work with PathCam imaging systems, Ensemble accessories, and formalin dispensing
- Frameless door
- Friction hinges

### Materials

- 304 stainless steel sheetmetal and hardware
- Lexan 3/8" clear polycarbonate plastic

### Dimensions

- 19.3" D x 17.6" H
- AX124L/R-048: 48.3" W
- AX124L/R-060: 60.3" W
- AX124L/R-072: 72.3" W

### Performance

- Ideal for open formalin processing operations
- Ensures maximum fume capture even in high traffic cross flow laboratory conditions

### Compatibility

- Fits the Mopec Maestro Encore grossing stations of the accompanying sized unit
- **Not** compatible with 96" units

### NOTE:

- This item **must** be ordered at the time of purchase to be properly installed during manufacturing of the grossing station.

## Variations

Model No.	Description	Dimensions
AX124L-048	Maestro Encore Safety Sash, Left Hand, 48" Units	48.3" W x 19.3" D x 17.6" H
AX124L-060	Maestro Encore Safety Sash, Left Hand, 60" Units	60.3" W x 19.3" D x 17.6" H
AX124L-072	Maestro Encore Safety Sash, Left Hand, 72" Units	72.3" W x 19.3" D x 17.6" H
AX124R-048	Maestro Encore Safety Sash, Right Hand, 48" Units	48.3" W x 19.3" D x 17.6" H
AX124R-060	Maestro Encore Safety Sash, Right Hand, 60" Units	60.3" W x 19.3" D x 17.6" H
AX124R-072	Maestro Encore Safety Sash, Right Hand, 72" Units	72.3" W x 19.3" D x 17.6" H

**MOPEC**  
GROUP

**mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI