

PathCam Micron Maestro Mount



SPECIFICATIONS

Construction

- Track-mountable articulating arm for vertical and horizontal adjustability designed to hold PathCam Micron Camera
- Threaded knobs vertically adjust arm mount along entire mounted track length

Materials

- Manufactured from a combination of resilient stainless steel and lightweight recycled aluminum

Capacity

- 5 lbs

Dimensions

- Articulating Camera arm:
 - 16.1" H x 17.9" – 22.3" L (depending on arm extension)
 - 180° of sweep and 13" of lift
- Focusing wheel adjusts camera mount vertically 10" on post
- Camera mount rotates 360° around post with locking set screw
- Ball and socket type peg allows for up to 45° of adjustable lateral viewing angles
- Rigid vertical peg provides a dedicated straight-on view

Performance

- Precision machine crafted components provide smooth adjustment and positioning
- Easy to use—simple installation and effortless operation
- Integrated cable management
- Balanced and stable ergonomic design

Compatibility

- Maestro and Maestro Encore grossing stations
- AX250 PathCam Camera Wall-Mount

NOTE:

- PathCam Micron Camera and PathCam Software sold separately

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI