

# Maestro Side Shields



## SPECIFICATIONS

### Construction

- Designed to mount on the outside edges of a Maestro or Maestro Encore grossing station
- Built to block outside disturbances, channel fumes into the grossing station vent grille, and act as splash guards
- Can be easily removed or replaced

### Materials

- Brackets: 304 stainless steel construction with #4 finish
- Hardware: 304 stainless steel
- Shields: clear polycarbonate

### Dimensions

- 20" W x 36.2" H x .25" Thick

### Compatibility

- Compatible with all Maestro and Maestro Encore grossing stations

**MOPEC**  
GROUP

**Mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI