

Maestro Dry Contact for Building Management System



SPECIFICATIONS

Construction

- Electrical wires to connect a Maestro or Maestro Encore grossing station to the facility Building Management System
- Signal wires run from facility to the grossing station and connect to the alarm relay signal output on electrical cabinet DIN rail
- Includes one set of dry contacts

Performance

- This feature connects to the facility's Building Management System to inform the appropriate parties that the grossing station has an alert condition
- BMS output triggering is setup in the Notifications Setting Page in the touchscreen panel

Compatibility

- Compatible with all Maestro and Maestro Encore grossing stations
- Maestro Encore units have this option integrated into the wiring harness and labeled "DRY CONTACTS DC1/DC2)"

NOTE:

- This option **must** be ordered at the time of purchase and installed during manufacturing of the grossing station.

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI