

Maestro 1/2 HP Disposal



SPECIFICATIONS

Construction

- Upgrade option for legacy grossing stations that were ordered without a disposal
- Now standard in all Maestro Encore grossing stations

Materials

- Polished stainless steel

Dimensions

- 7.5" L x 7.5" W x 13.375" H

Performance

- 1/2 Horsepower
- 2600 RPM

Compatibility

- Compatible with all Maestro and Maestro Encore grossing stations

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI