

Maestro 2X Quick Coupler Base Magnifier



SPECIFICATIONS

Construction

- Designed to mount to the shelf of your Maestro or Maestro Encore grossing station to provide 2X magnification
- Includes a coupler base for quick-release swivel mounting

Materials

- Glass 2X magnifier
- Vinyl-coated flex arm

Dimensions

- Magnifier: 4" diameter
- Flexible Arm: 24" L

Compatibility

- Compatible with all Maestro and Maestro Encore grossing stations

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI