

Maestro Non-Elevating Leg Frame



MB1030-048



MB1030-060



MB1030-072



MB1030-096

SPECIFICATIONS

Construction

- Designed for customers who do not need the elevation function on a standard Elevating Maestro or Maestro Encore grossing station
- Includes leg frame and leveling feet with 1.625" of adjustment
- If you'd like to add to an existing Countertop Maestro or Maestro Encore grossing station, please contact your Mopec sales representative

Materials

- 304 Stainless steel, 180 grit finish
- 11 ga 304 stainless steel seismic anchor brackets included, with 3/8" 304 stainless hardware

Dimensions

- 29" D x 32" H
- MB1030-048: 48" L
- MB1030-060: 60" L
- MB1030-072: 72" L
- MB1030-096: 96" L

Compatibility

- Compatible with all Maestro and Maestro Encore Non-Elevating grossing stations
- MB1030-048: 48" stations
- MB1030-060: 60" stations
- MB1030-072: 72" stations
- MB1030-096: 96" stations

Variations

Model No.	Description	Dimensions
MB1023-048	Non-Elevating Leg Frame, 48" Units	48" L x 29" D x 32" H
MB1023-060	Non-Elevating Leg Frame, 60" Units	60" L x 29" D x 32" H
MB1023-072	Non-Elevating Leg Frame, 72" Units	72" L x 29" D x 32" H
MB1023-096	Non-Elevating Leg Frame, 96" Units	96" L x 29" D x 32" H

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI