



# *Grossing Consumables Toolkits*



# Pathologists' Assistant Approved Grossing Consumables Toolkits

*Mopec's grossing toolkits have been created with the expertise and experience of Pathologists' Assistants to address the basic needs in every gross room. Simplify your grossing consumables management with a subscription to Mopec's basic toolkits. All consumables are Mopec Certified for quality, performance, and compatibility.*

## FORMALIN MANAGEMENT TOOLKIT #BU005

Formalin and formalin neutralizing in one simple package. This toolkit includes the standard 2.5 gallon cube of formalin along with neutralizer pads for safer cutting and neutralizer powder for safer disposal of spent formalin. *Mopec Formalin is a 10% neutral buffered formalin used as a general purpose fixative for tissues. While Mopec formalin penetrates quickly, it fixes slowly. No other dilution or addition necessary. Mopec Formalin is also available in prefilled specimen containers or large, bulk containers.*

PRODUCT #	PRODUCT DESCRIPTION	QTY	FREQUENCY
BE094-CS	FORMAWAY NEUTRALIZING PADS, 8" x 11", 100 PACK	2	2 Months
BG550-12	FORMAWAY NEUTRALIZER POWDER, CASE OF 12	1	2 Months
BG520	MOPEC PREFILLED FORMALIN CUBE, 2.5 GALLONS	4	2 Months



## CLEANING & MAINTENANCE TOOLKIT #BU004

Keep your grossing station like new with Mopec's cleaning & maintenance toolkit. This toolkit comes complete with all of the items you need to properly care for your Mopec Maestro® Encore and your Mopec dissection instruments. *Consult the Mopec Maestro user manual for proper cleaning and maintenance procedures.*



PRODUCT #	PRODUCT DESCRIPTION	QTY	FREQUENCY
BE045	SANIPATH DISINFECTANT FOAMING SPRAY	1	3 Months
BE046	EAZYSOAK ENZYME INSTRUMENT CLEANER	1	3 Months
BE048	STAINLESS STEEL CLEANER AND POLISH	1	3 Months
BE042	PATHODORIZER+, ODOR ELIMINATING SOLUTION	1	3 Months
BF046	MAESTRO STANDARD FILTER, FOR 60" MAESTROS	1	3 Months

*Other variations available for 48", 72", and 96" units*

## LARGE DISSECTION TOOLKIT #BU003

Designed for Pathologists' Assistants primarily performing dissections on large specimens, the large dissection toolkit includes the essentials for safe grossing with a #60 scalpel blade. Stop removing scalpel blades with fingers or forceps with the included Bladeflask EVO Blade Remover. Extra large Formaway Formalin Neutralizing Pads ensure your formalin exposure is reduced while also keeping your scalpel blades sharper longer.

PRODUCT #	PRODUCT DESCRIPTION	QTY	FREQUENCY
AJ112	SCALPEL BLADES, STERILE, #60 - 100 PACK	2	Monthly
AJ160	BLADEFLASK EVO UNIVERSAL BLADE REMOVER	2	Monthly
BE094	FORMAWAY NEUTRALIZING PADS, 11" x 15", 25 PACK	2	Monthly





## BIOPSY & SMALL DISSECTION TOOLKIT #BU002

Designed for Pathologists' Assistants and Grossing Technicians primarily performing routine biopsies and dissections on small specimens, this toolkit includes the essentials for safe grossing with a #22 scalpel blade. In addition to the BladeFlask EVO Blade Remover and small Formaway Formalin Neutralizing Pads, specimen protection and loss prevention is covered with Mopec Premium Biopsy Bags, biopsy wraps, and papers.

PRODUCT #	PRODUCT DESCRIPTION	QTY	FREQUENCY
AJ105	SCALPEL BLADES, STERILE, #22 - 100 PACK	2	Monthly
AJ160	BLADEFLASK EVO UNIVERSAL BLADE REMOVER	2	Monthly
BE093	FORMAWAY NEUTRALIZING PADS, 8" x 11", 25 PACK	2	Monthly
BI004	PREMIUM BIOPSY BAGS, SMALL, 1000/BOX	1	Monthly
BI056	BxPADS, BLUE SPONGE - 100 PACK	1	Monthly
BI059	BxWRAPS, 500 PACK, 2" x 3"	1	Monthly



## ALL-IN-ONE TOOLKIT #BU001

The All-in-One Toolkit delivers everything Pathologists' Assistants and Grossing Technicians need for safe, efficient, and compliant grossing. This versatile kit combines the key elements from our small specimen, large specimen, cleaning & maintenance, and formalin management toolkits to support every aspect of your workflow.

PRODUCT #	PRODUCT DESCRIPTION	QTY	FREQUENCY
AJ105	SCALPEL BLADES, STERILE #22, 100 PACK	4	3 Months
AJ112	SCALPEL BLADES, STERILE #60, 100 PACK	4	3 Months
AJ160	BLADEFLASK EVO UNIVERSAL BLADE REMOVER	4	3 Months
BI004	PREMIUM BIOPSY BAGS, SMALL, 1000 PACK	3	3 Months
BI056	BXPADS, BLUE SPONGE, 100 PACK	2	3 Months
BE093-CS	FORMAWAY FORMALIN. NEUTRALIZING PADS, 8" X 11", 100 PACK	2	3 Months
BE094-CS	FORMAWAY FORMALIN. NEUTRALIZING PADS, 11" X 15", 100 PACK	2	3 Months
BE045	SANIPATH DISINFECTANT FOAMING SPRAY	1	3 Months
BE046	EAZYSOAK INSTRUMENT CLEANER	1	3 Months
BE048	CLEARSTEEL STAINLESS STEEL CLEANER & POLISH SPRAY	1	3 Months
BE042	PATHODORIZER PLUS, DEODORIZER, CLEANER & DEGREASER	1	3 Months
BF046	MAESTRO STANDARD FILTER, FOR 60" MAESTROS	1	3 Months

Other variations available for 48", 72", and 96" units



## Simplify your Grossing Consumables Management!

SCAN QR CODE OR VISIT [MOPEC.COM/TOOLKITS](https://mopec.com/toolkits)

# Other *Mopec Certified* Grossing Consumables

*Mopec's comprehensive range of mission-critical consumables are designed to support every stage of anatomic pathology, from grossing to diagnosis. As a leading manufacturer and distributor, Mopec delivers top-quality products essential for your laboratory needs.*


## *Tissue Cassettes*

Designed to meet the needs and requirements of our customers, our line of Mopec Cassettes are compatible with a wide range of tissue cassette printers. Use our Tissue Cassette Selection Tool to choose your preferred lid style, cassette type, packaging, and color to fit your exact tissue printer.



**Free samples available!**

SCAN QR CODE OR VISIT [MOPEC.COM/CASSETTES](https://mopec.com/cassettes)



## *Harmony*

TISSUE MARKING DYE

Mopec's Harmony dyes are used to permanently mark surgical tissue margins in Histological specimens. Marking dyes are available in 9 colors, with a standardized viscosity over the range. This ensures that there is less variation in viscosity from one color to another as compared to competitor products and allows the easy marking of both small and large specimens with any color.



**Free samples available!**

SCAN QR CODE OR VISIT [MOPEC.COM/HARMONY](https://mopec.com/harmony)



## *Scalpels & Knives*

Reduce your risk to injury and infectious diseases with Mopec scalpels and knives. This line includes scalpel handles, scalpel blades, and knives in various sizes and materials. Mopec scalpels and knives are designed to keep your employees safe and fit the needs of your lab.



SCAN QR CODE OR VISIT [MOPEC.COM/SCALPELS](https://mopec.com/scalpels)



Shop All Anatomic Pathology Consumables And Supplies

SCAN QR CODE OR VISIT [MOPEC.COM/SUPPLIES](https://mopec.com/supplies)