


# Maestro Red GFCI and USB Outlets



## SPECIFICATIONS

### Construction

- Designed for customers who need red outlets to indicate they are connected to a backup power source
- Includes (1) heavy-duty, hospital-grade, tamper-resistant, duplex receptacle outlet with two type-A USB ports; (1) extra heavy-duty, hospital grade, self-test GFCI duplex receptacle outlet; and wall plate

### Performance

- Features power indication, 15 amp and 20 amp feed-through, and 125 volt
- Receptacle outlet with USB is located inside the Maestro column behind touchscreen
- GFCI receptacle outlet is located on the outside of the Maestro column on the opposite side of sink

### Compatibility

- Compatible with all Maestro and Maestro Encore grossing stations
- Compatible with back or side wiring, NEMA 5-15R, 2-pole, 3-wire, and self-grounding

**MOPEC**  
GROUP

**Mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI