

Maestro Formalin Dispensing Tap for FP100 Formalin Pumping Station



SPECIFICATIONS

Construction

- Designed to be used on Maestro and Maestro Encore grossing stations in conjunction with a centrally-located FP100 Formalin Pumping Station
- Includes (1) dispensing tap

Materials

- Chrome-plated brass

Compatibility

- Compatible with all Maestro and Maestro Encore grossing stations
- To be used in conjunction with the FP100 Formalin Pumping Station

NOTE:

- This option is different than the MB1040 Formalin Dispensing System and the MB1043 Formalin Dispensing and Collection System, which have an internal tank within the grossing station for pumping formalin.
- This must option **must** be used with an external formalin pumping station, such as the FP100

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI