

Maestro Formalin Dispensing



SPECIFICATIONS

Construction

- Designed to be used on Maestro and Maestro Encore grossing stations to dispense formalin from an internal dispensing tank
- Includes (1) dispensing tap with flow control, formalin transfer pump, and internal dispensing tank

Performance

- Gravity-fed dispensing system utilizes a siphon hose and transfer pump to fill a reservoir tank located in the hood of the unit from facility-provided formalin supply container
- Equipped with level sensors that turn off the pump when full
- Alarm will alert users when tank becomes full, or when tank becomes low

Materials

- Spigot: Chrome-plated brass

Capacity

- 5 gallons

Compatibility

- Compatible with all Maestro and Maestro Encore grossing stations

NOTE:

- This item **must** be ordered at the time of purchase to be properly installed during manufacturing of the grossing station.
- This option is different than the MB1037 Formalin Dispensing Tap, which must be used with an external formalin pumping station, such as the FP100
- Formalin **not** included

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI