

Maestro Formalin Collection



SPECIFICATIONS

Construction

- Designed to be used on Maestro and Maestro Encore grossing stations to collect formalin into a collection carboy
- Includes (1) formalin collection funnel, nylon hose, collection carboy with level sensor, and collection caddy with leak sensor

Performance

- Gravity-fed collection system that drains from the funnel in the sink into the carboy situated in the caddy
- Carboy equipped with a sensor and alarm to indicate when system is full

Capacity

- 2.5 gallons

Compatibility

- Compatible with all Maestro and Maestro Encore grossing stations
- BA077: Funnel with attached screen

NOTE:

- Formalin **not** included

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI