


Maestro Infrared Faucet



SPECIFICATIONS

Construction

- Designed to control the On/Off function for the sink on a Maestro or Maestro Encore grossing station
- Includes infrared sensor, water solenoid, and coil

Performance

- Controls the On/Off function of the unit's sink
- Wave hand within sensor area to turn on faucet
- Faucet will shut off automatically if movement is no longer detected, or after a designated time set through water timer

Compatibility

- Compatible with all Maestro and Maestro Encore grossing stations

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI