



REV	DESCRIPTION	DATE	ECN #	APPR
A	INITIAL RELEASE	00-XXX-19	XXXXX	MM

Item Number	Quantity	File Name	Title	Material
1	1	PTE-0001-26D.asm	PTA-45 ARM	
2	1	S7972.asm	PIVOT ADAPTER / CABLE PASS THROUGH	
3	1	AX250HD-001.asm	HD CAMERA ENCLOSURE	
4	1	91280A351_Medium-Strength Class 8.8 Steel Hex Head Screw.par		
5	1	6092K21_Nylon Seven Arm Knob.par		
6	4	90600A224_18-8 Stainless Steel Button Head Torx Screws.par	SS Button Head Torx Screw, 10-32 X 5/16"	
7	1	PM0766.par	ERGOTRON 42" WALL TRACK	
8	2	S6836-B.par	AX ACCESSORY MOUNT	RM0258, FLAT, AL, 1/2" X 2-1/2", 6061
9	1	PM0354-B.asm	ERGOTRON, KEYBOARD ARM	
10	1	PM0355-B.asm	ERGOTRON, MONITOR ARM	
11	1	PM0875.par	ERGOTRON 10" WALL TRACK	
12	1	PM0864.par	ERGOTRON, WALL TRACK MOUNT BRACKET KIT FOR CPU HOLDER	
13	1	PM0874.asm	ERGOTRON, CPU HOLDER	
14	4	PL0091.par	5/16"-18 x 1-1/2"L KNOB	Stainless steel

THIS DOCUMENT CONTAINS PROPRIETARY INFORMATION OF MOPEC, AND MAY NOT BE COPIED IN WHOLE OR PART, AND NO USE SHALL BE MADE OF IT OTHER THAN FOR WHICH IT HAS BEEN SUBMITTED WITHOUT THE WRITTEN PERMISSION OF A MOPEC MANAGER

MOPEC 2021

UNLESS OTHERWISE SPECIFIED, ALL DIMENSIONS IN INCHES

TOLERANCES

0.0	± 0.1 in
0.00	± 0.05 in
0.000	± 0.010 in
0.0000	± 0.0010 in
ANGLES	±0.5°

BREAK ALL CORNERS 0.063 X 45°

ALL FILLET RADI R0.063 /0.125in

INSPECT OVAL DIMENSIONS

MATERIAL

**MOPEC GROUP**  
Advanced Technology Manufacturing

**MORTECH**  
MANUFACTURING

DRAWN TO 3RD ANGLE PROJECTION

DRAWN BY: dstewart2  
DATE: 5/29/2025

CHECKED BY:   
DATE: 5/29/2025

APPROVED BY: MMcCLAIN  
DATE: 5/29/2025

PART NUMBER: AX250HD-000	
DESCRIPTION: PathCam Lite	
SIZE: A2	FILE NAME: AX250HD-000.dft
SCALE: 1:4	DO NOT SCALE IF IN DOUBT ASK
SHEET 1 OF 2	



