

PathCam Wall Mount System

The PathCam wall mount system enhances pathology workflow by integrating an articulating keyboard arm, monitor arm, and camera arm for ergonomic, space-saving digital documentation. Designed for use at grossing stations, it improves efficiency, accuracy, and collaboration in clinical settings through seamless integration of digital systems.



AX250HD



AX250HR

SPECIFICATIONS

Construction

- (1) 42" Wall-mount track
- (1) 10" Wall-mount track
- (1) Articulating Camera arm with 60" reach
- (1) Vertical CPU holder for 10" wall track
- Mounting hardware for wall tracks included

Materials

- Powder-coated cast aluminum grade 6061
- 304 stainless steel construction with a polished No. 4 finish

Dimensions

- Articulating Camera Arm:
 - Maximum Reach: 62" (wall surface to camera-mounting mid-point)
 - Minimum Reach: 20"
 - Rear pivot to mid-point of arm: 31"
 - Vertical CPU holder size Range: 1-3/8" to 3-3/4"
- Keyboard Arm: 17.9" – 22.3" L x 33.5" W x 16.1" H
- Monitor Arm: 12.3" L x 10.2" W x 17.2" H

Performance

- Has several pivot points to articulate up and down, and side to side.
- Can pivot at the wall, swivel, and tilt to optimize camera-viewing angle
- Can be pulled out when in use, then pushed back when not in use.

Compatibility

- Compatible with PathCam HD and HR Cameras
- AX032 Keyboard Arm
- AX034 Monitor Arm
- AX047 PathCam Foot Pedal Command

Variations

| Model No. | Description | Compatibility |
|-----------|------------------------------|-------------------|
| AX250HR | PathCam HR Wall Mount System | PathCam HR Camera |
| AX250HD | PathCam HD Wall Mount System | PathCam HD Camera |

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI