

PathCamHD Retrofit (Camera Only)

The Mopec AX265 PathCamHD Retrofit (Camera Only) enhances any Maestro Encore grossing stations with high-definition imaging and streaming capabilities, ideal for teaching and demonstrations. Designed for intuitive use, it comes preprogrammed with PathCam software settings for consistent visuals. Its compact setup ensures minimal disruption to workflow and a seamless solution for pathology environments.



SPECIFICATIONS

Construction

- Designed for those who are using PathCam Systems for live demonstrations or teaching/streaming applications on their legacy grossing station
- Preprogrammed with PathCam Software settings
- Includes PathCam HD Camera, HDMI-USB hub converter, HDMI cable, 60-day PathCam Software free trial, and 1-year limited manufacturer's warranty

Performance

- Image Sensor: 2.1 MP CMOS
- Live video feed: 1920 x 1080 pixels
- Captures: 0.2–4.5 MB PNG files; Full HD 146–202 MB/min AVI files
- Optical/digital magnification 10x/12x
- Max FOV: 11.3 x 8 ft (3.4 x 2.4 m) @ 9 ft (2.7 m)
- Min FOV: 7.5 x 5.3 in (19.1 x 13.5 cm) @ 6 in (15.2 cm)
- Auto Focus
- Auto White Balance
- Auto Exposure
- Multi-monitor switching system compatible

Compatibility

- Mopec Maestro grossing stations
- Legacy (MB, ML, ME, MG) grossing stations
- See 'NOTE' section for details on correct Retrofit Installation Kit required for models, not included in purchase of AX265

NOTE:

- Requires a Retrofit Installation Kit (AX043) based on existing equipment
- 2018–2022 Maestro Platform: AX043-L1
- 2023–2024 Maestro Platform: AX043-L2
- Legacy Stations (MB, ML, ME, MG): AX043-L3
- PathCam Software sold separately

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI