


PathCamHR Retrofit (Camera Only)

The Mopec AX266 PathCamHR Retrofit (Camera Only) enhances any Maestro Encore grossing stations with high-definition imaging and streaming capabilities, ideal for teaching and demonstrations. Designed for intuitive use, it comes preprogrammed with PathCam software settings for consistent visuals. Its compact setup ensures minimal disruption to workflow and a seamless solution for pathology environments.



SPECIFICATIONS

Construction

- Designed for those who are using PathCam Systems for diagnostic recording and publishing on their legacy grossing station
- Preprogrammed with PathCam Software settings
- Includes PathCam HR 6K Camera, power USB cable kit, camera SD card, 60-day PathCam Software free trial, and 1-year limited manufacturer's warranty

Performance

- Image Sensor: 24 MP CMOS
- Live video feed: 640 x 480 pixels
- Captures: 0.1–60 MB PNG files; Full HD 340 MB/min MOV files
- Optical/digital magnification: 10x/15x
- Max FOV: 45 x 32 in (114.3 x 81.3 cm) @ 3 ft (91.4 cm)
- Auto Focus
- Auto White Balance
- Auto Exposure

Compatibility

- Mopec Maestro grossing stations
- Legacy (MB, ML, ME, MG) grossing stations
- See 'NOTE' section for details on correct Retrofit Installation Kit required for models, not included in purchase of AX266

NOTE:

- Requires a Retrofit Installation Kit (AX043) based on existing equipment
- 2018–2022 Maestro Platform: AX043-L1
- 2023–2024 Maestro Platform: AX043-L2
- Legacy Stations (MB, ML, ME, MG): AX043-L3
- PathCam Software sold separately

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI