


Dissection Cart with Fluid Collection and Removable 27" Body Tray

Mopec's Dissection Cart with Fluid Collection offers a streamlined and hygienic solution for pathology and anatomy labs. Elevate your workflow with this ergonomic and durable solution, crafted to enhance safety, cleanliness, and productivity in your laboratory environment.



SPECIFICATIONS

Construction

- Stainless steel construction
- (4) heavy duty 5" casters with (3) total lock and (1) directional lock
- (4) Post sockets included
- Conveyor roller-transfer style top
- Built-in slope
- Lower shelf
- 2.5-gallon carboy container and carrier

Materials

- 304 stainless steel construction with a polished No. 4 finish
- Welded stainless steel type 304 11GA frame

Dimensions

- 80" L x 30.9" W x 35.3" – 35.8" H

Capacity

- 350 Lbs

Options & Accessories

- JC027: 27" tray without drain hole (sold separately)
- JC027D: 27" tray with drain hole, rubber plug, and stainless steel sash chain (sold separately)
- HO005: Post leg-holders
- HO006: Book holder
- HO030: Instrument tray holder

NOTE:

- Tray **not** included. JC027 and JC027D trays sold separately.

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI