

Pathport 3 Specimen Transport & Storage

The PathPort 3 Specimen Transport & Storage System is the ultimate solution for safely transporting and temporarily storing histology and pathology specimens. Designed with durability, security, and efficiency in mind, the PathPort 3 was developed to address the critical need for a secure, compliant, and efficient specimen transport system in hospital and laboratory settings.



SPECIFICATIONS

Construction

- Heliarc-welded seams and joints, ground and polished to match adjoining surfaces
- Includes safety locks and removable bottom shelf

Materials

- 304 stainless steel construction with a polished No. 4 finish

Dimensions

- Overall: 40.1" L x 18.1" W x 42.4" H
- Inside: 34" L x 16" W x 11.6" H

Performance

- OSHA Compliant
- Equipped with a Pathport fitted formaldehyde neutralizing pad

Capacity

- 80 lbs

Options & Accessories

- BE093: Replacement Disposable Pads (25 pack)
- Note: This will require stacking 2, side by side, in your Pathport

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI